

Product Datasheet

Malaria (pf-HRP-2) Antibody HM245

Product Name	Malaria (pf-HRP-2) Antibody
Catalog Number	HM245
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Note	Detection capability of ≥ 50 parasites/ml
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 %
Buffer	20 mM PBS, 50 mM NaCl, pH 7.4
Cross Reactivity	No reactivity with other Plasmodium species
Specificity	Plasmodium falciparum (malaria) HRP-2
Matched Pairs	HM246 (capture)
Preservatives	0.1% Sodium Azide