

Product Datasheet

Malaria (pf-HRP-2) Antibody HM246

Product Name	Malaria (pf-HRP-2) Antibody
Catalog Number	HM246
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Detection capability of ≥ 50 parasites/ml
Isotype	IgG1
Clonality	Monoclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 8.0 with 50 mM NaCl
Cross Reactivity	No reactivity with other Plasmodium species
Specificity	Plasmodium falciparum (malaria) HRP-2
Matched Pairs	HM245 (detection)
Preservatives	0.1% Sodium Azide