

Product Datasheet

Malaria (pan-pLDH) Antibody HM249

Product Name	Malaria (pan-pLDH) Antibody
Catalog Number	HM249
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Recombinant Plasmodium LDH
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Note	Detection capability of ≥ 50 parasites/ml
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 %
Buffer	10 mM Phosphate Buffered Saline, pH 8.0 with 50 mM NaCl
Specificity	Malaria pan-pLDH: P. falciparum, P. vivax, P. ovale, P. malariae
Matched Pairs	HM250 (capture)
Preservatives	0.1% Sodium Azide