

Product Datasheet

Malaria (pan-pLDH) Antibody HM250

Product Name	Malaria (pan-pLDH) Antibody
Catalog Number	HM250
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Recombinant Plasmodium LDH
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Note	Detection capability of ≥ 50 parasites/ml
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 %
Buffer	20 mM Phosphate Buffered Saline, pH 7.4 with 50 mM NaCl
Specificity	Malaria pLDH
Matched Pairs	HM249 (detection)
Preservatives	0.1% Sodium Azide