

Product Datasheet

Hepatitis B Surface Antigen (HBsAg) Antibody HM252

Product Name	Hepatitis B Surface Antigen (HBsAg) Antibody
Catalog Number	HM252
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Note	Sensitivity of 0.5 ng/ml in Lateral Flow
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 %
Buffer	10 mM Phosphate Buffered Saline, pH 8.0 with 50 mM NaCl
Specificity	HBsAg (including adw, ayw, and adr subtypes)
Matched Pairs	HM251 (detection)
Preservatives	0.1% Sodium Azide