

Product Datasheet

Hepatitis B Surface Antigen (HBsAg) Antibody HM252

Product Name	Hepatitis B Surface Antigen (HBsAg) Antibody
Catalog Number	HM252
Description	Ion Exchange Purified Monoclonal Antibody
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Note	Sensitivity of 0.5 ng/ml in Lateral Flow
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	≥ 95 %
Buffer	20 mM phosphate-buffered saline, pH 7.4 ± 0.2, with 50 mM NaCl, 0.1% sodium azide as a preservative.
Specificity	Hepatitis B surface antigen
Matched Pairs	HM251 (detection)