

Product Datasheet

Neisseria gonorrhoeae Antibody HM253

Product Name	Neisseria gonorrhoeae Antibody
Catalog Number	HM253
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Note	Detection capability of $\leq 10^4$ cfu/ml
Isotype	IgG2b
Clonality	Monoclonal
Purity	> 95 %
Buffer	10 mM Phosphate Buffered Saline, pH 8.0 with 50 mM NaCl
Cross Reactivity	No reactivity with <i>N. meningitidis</i> , <i>N. cinerea</i> , <i>N. lactamica</i> , <i>M. sicca</i> , <i>B. catarrhalis</i> , <i>E. coli</i> , <i>P. mirabilis</i> , <i>G. vaginalis</i> , <i>C. trachomatis</i>
Specificity	<i>N. gonorrhoeae</i>
Matched Pairs	HM254 (capture)
Preservatives	0.1% Sodium Azide