

Product Datasheet

Hemoglobin (Hb) Antibody HM257

Product Name	Hemoglobin (Hb) Antibody
Catalog Number	HM257
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Note	Sensitivity of ≤ 30 ng/ml
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 %
Buffer	10 mM Phosphate Buffered Saline, pH 8.0 with 50 mM NaCl
Cross Reactivity	No reactivity with hemoglobin from pig, cow, goat, sheep, horse, chicken, rabbit
Specificity	Human hemoglobin
Matched Pairs	HM258 (capture)
Preservatives	0.1% Sodium Azide