

## Product Datasheet

### Aflatoxin B1 (AFT B1) Antibody HM260

<b>Product Name</b>	Aflatoxin B1 (AFT B1) Antibody
<b>Catalog Number</b>	HM260
<b>Description</b>	Protein A Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Monoclonal Antibody
<b>Storage</b>	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Note</b>	Meets 5 ng/ml cut off
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95 %
<b>Buffer</b>	20 mM Phosphate Buffered Saline, pH 7.4 with 50 mM NaCl
<b>Cross Reactivity</b>	No reactivity with Aflatoxins M1 or M2
<b>Specificity</b>	Aflatoxin B1, G1, Q1, B2, G2
<b>Matched Pairs</b>	LA441 (capture)
<b>Preservatives</b>	0.1% Sodium Azide