

Product Datasheet

Aflatoxin B1 (AFT B1) Antibody HM260

| | |
|----------------------------|--|
| Product Name | Aflatoxin B1 (AFT B1) Antibody |
| Catalog Number | HM260 |
| Description | Protein A Purified Monoclonal Antibody |
| Species/Host | Mouse |
| Tested Applications | ELISA, LF |
| Brand | EastCoast Bio |
| Form/Appearance | Protein A Purified Monoclonal Antibody |
| Storage | Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| Note | Meets 5 ng/ml cut off |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | > 95 % |
| Buffer | 20 mM Phosphate Buffered Saline, pH 7.4 with 50 mM NaCl |
| Cross Reactivity | No reactivity with Aflatoxins M1 or M2 |
| Specificity | Aflatoxin B1, G1, Q1, B2, G2 |
| Matched Pairs | LA441 (capture) |
| Preservatives | 0.1% Sodium Azide |