

Product Datasheet

Estradiol (E2) Antibody HM275

Product Name	Estradiol (E2) Antibody
Catalog Number	HM275
Description	Protein G Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Protein G Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 85 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	Estradiol
Preservatives	0.03% Proclin 300