

## Product Datasheet

### Estradiol (E2) Antibody HM275

<b>Product Name</b>	Estradiol (E2) Antibody
<b>Catalog Number</b>	HM275
<b>Description</b>	Protein G Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein G Purified Monoclonal Antibody
<b>Storage</b>	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 85 % by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Specificity</b>	Estradiol
<b>Preservatives</b>	0.03% Proclin 300