

## Product Datasheet

### Methcathinone Antibody HM300

<b>Product Name</b>	Methcathinone Antibody
<b>Catalog Number</b>	HM300
<b>Description</b>	Protein G Purified Monoclonal Antibody
<b>Immunogen</b>	Methcathinone-BSA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein G purified monoclonal antibody
<b>Storage</b>	Short term: 2-8 °C. Long term: -20 °C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	≥ 95 % by SDS-PAGE and HPLC
<b>Buffer</b>	0.015 M potassium phosphate, pH 7.2 ± 0.2, with 0.85% sodium chloride, 0.05% sodium azide as a preservative.
<b>Specificity</b>	Methcathinone