

Product Datasheet

Ferritin Antibody HM302

| | |
|------------------------|--|
| Product Name | Ferritin Antibody |
| Catalog Number | HM302 |
| Description | DEAE Purified Monoclonal Antibody |
| Immunogen | Ferritin isolated from human liver |
| Brand | EastCoast Bio |
| Form/Appearance | DEAE Purified Monoclonal Antibody |
| Storage | Short Term: 2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| Isotype | Mouse IgG1 |
| Clonality | Monoclonal |
| Purity | ≥ 95 % by SDS-PAGE and HPLC |
| Buffer | 0.015 M potassium phosphate, pH 7.2 ± 0.2, with 0.85% sodium chloride, 0.05% sodium azide as a preservative. |
| Specificity | Human Ferritin |
| Matched Pairs | HM303 (detection) |