

Product Datasheet

Ferritin Antibody HM303

Product Name	Ferritin Antibody
Catalog Number	HM303
Description	DEAE Purified Monoclonal Antibody
Immunogen	Ferritin isolated from human liver
Brand	EastCoast Bio
Form/Appearance	DEAE purified monoclonal antibody
Storage	Short Term: 2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	≥ 95 % by SDS-PAGE and HPLC
Buffer	0.015 M potassium phosphate, pH 7.2 ± 0.2, with 0.85% sodium chloride, 0.05% sodium azide as a preservative.
Specificity	Human ferritin
Matched Pairs	HM302 or HM304 (capture or detection)