

Product Datasheet

Rabies Virus Antibody HM305

Product Name	Rabies Virus Antibody
Catalog Number	HM305
Species/Host	Mouse
Tested Applications	ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline pH 7.2
Cross Reactivity	Reactive with many Rabies Virus Strains
Specificity	Rabies Virus surface glycoprotein
Preservatives	0.1% Sodium Azide