

Product Datasheet

Neutrophil Gelatinase-associated Lipocalin (NGAL) Antibody HM307

Product Name	Neutrophil Gelatinase-associated Lipocalin (NGAL) Antibody
Catalog Number	HM307
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Purified native human NGAL
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): $+2-8$ °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	≥ 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	NGAL
Matched Pairs	HM306 (capture)
Preservatives	0.09% Sodium Azide