

Product Datasheet

Zopiclone Antibody HM310

Product Name	Zopiclone Antibody
Catalog Number	HM310
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Zopiclone-BSA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE and HPLC
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	No reactivity with Zolpidem
Specificity	Zopiclone
Matched Pairs	LA269 (capture)
Preservatives	0.05% Sodium Azide