

Product Datasheet

Human CD141 ELISA Kit 850.720.192

Product Name	Human CD141 ELISA Kit
Catalog Number	850.720.192
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	0.6-20 ng/ml
Sensitivity	0.38 ng/ml
Tested Applications	ELISA
Alternative Names	Thrombomodulin
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human soluble molecules
Specificity	Recognizes both natural and recombinant human soluble thrombomodulin
Sample Types	Serum, Plasma, Cell culture supernatant