

Product Datasheet

Zolpidem Antibody HM312

Product Name	Zolpidem Antibody
Catalog Number	HM312
Description	DEAE Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Zolpidem-BSA
Brand	EastCoast Bio
Form/Appearance	DEAE Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE and HPLC
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	Zolpidem Urinary Antigen (Phenyl-4-Carboxylic Acid): 100%, Zolpidem: 90%
Specificity	Zolpidem
Matched Pairs	LA155 (capture)
Preservatives	0.05% Sodium Azide