

Product Datasheet

NT-proBNP (61-76) Antibody HM316

Product Name	NT-proBNP (61-76) Antibody
Catalog Number	HM316
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF, IFA
Immunogen	Purified Recombinant Human NT-proBNP
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	NT-proBNP
Matched Pairs	HM315 (capture)
Preservatives	0.09% Sodium Azide