

Product Datasheet

Calprotectin Antibody HM320

Product Name	Calprotectin Antibody
Catalog Number	HM320
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Purified Human calprotectin
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	Phosphate Buffered Saline, pH 7.0
Specificity	Native Human Calprotectin
Matched Pairs	HM319 (capture)
Preservatives	0.09% Sodium Azide