

Product Datasheet

Pepsinogen I (PGI) Antibody HM322

Product Name	Pepsinogen I (PGI) Antibody
Catalog Number	HM322
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Purified Human pepsinogen I
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	PGI
Matched Pairs	HM321 (capture)
Preservatives	0.05% Sodium Azide