

Product Datasheet

Human CD21 ELISA Kit 850.730.192

Product Name	Human CD21 ELISA Kit
Catalog Number	850.730.192
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	3.12-100 U/ml
Sensitivity	3 U/ml
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human soluble molecules
Specificity	Recognizes human soluble CD21
Sample Types	Serum, Plasma, Cell culture supernatant