

Product Datasheet

Tramadol Antibody HM342

Product Name	Tramadol Antibody
Catalog Number	HM342
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Tramadol-KLH
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Note	Meets 100 ng/ml in urine testing
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Specificity	Tramadol
Preservatives	0.05% Sodium Azide