

Product Datasheet

Oxycodone Antibody HM343

Product Name	Oxycodone Antibody
Catalog Number	HM343
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Oxycodone-BSA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Note	Meets 100 ng/ml in urine testing
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	No reactivity with codeine, hydrocodone or hydromorphone and = 1% with morphine
Specificity	Oxycodone
Matched Pairs	LA266 (capture)
Preservatives	0.05% Sodium Azide