

Product Datasheet

Human CD95 ELISA Kit 850.740.192

Product Name	Human CD95 ELISA Kit
Catalog Number	850.740.192
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	93.75-3000 pg/ml
Sensitivity	47 pg/ml
Tested Applications	ELISA
Alternative Names	APO-1, FAS
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human soluble molecules
Specificity	Recognizes both natural and recombinant soluble human Fas
Sample Types	Serum, Plasma, Cell culture supernatant