

Product Datasheet

Dengue Virus (DENV-Env) Antibody HM348

Product Name	Dengue Virus (DENV-Env) Antibody
Catalog Number	HM348
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Recombinant Dengue Virus Envelope Protein
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	No reactivity with Zika or Chikungunya viruses
Specificity	All Dengue Virus serotypes
Preservatives	0.05% Sodium Azide