

Product Datasheet

Benzodiazepine (Oxazepam) Antibody HM350

Product Name	Benzodiazepine (Oxazepam) Antibody
Catalog Number	HM350
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Oxazepam-BSA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl
Cross Reactivity	Minimal reactivity with Clonazepam, Flurazepam, Prazepam
Specificity	Oxazepam and related compounds including Alprazolam, Diazepam, Lorazepam, Nitrazepam and Tremazepam
Matched Pairs	LA264 (capture)
Preservatives	0.05% Sodium Azide