

Product Datasheet

Human CD178 ELISA Kit 850.750.096

| | |
|----------------------------|--|
| Product Name | Human CD178 ELISA Kit |
| Catalog Number | 850.750.096 |
| Description | Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure. |
| Range | 62.5-2000 pg/ml |
| Sensitivity | 12 pg/ml |
| Tested Applications | ELISA |
| Alternative Names | FasL, CD95L |
| Brand | Diaclone |
| Storage | Store at +4°C |
| Cross Reactivity | No cross reactivity with other human soluble molecules |
| Specificity | Recognizes both natural and recombinant soluble human Fas ligand |
| Sample Types | Serum, Plasma, Cell culture supernatant |