

Product Datasheet

Procalcitonin (PCT) Antibody HM364

Product Name	Procalcitonin (PCT) Antibody
Catalog Number	HM364
Species/Host	Mouse
Tested Applications	LF, ELISA, WB
Immunogen	Recombinant Procalcitonin
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Procalcitonin
Matched Pairs	HM365 (capture)
Preservatives	0.05% Proclin 300