

Product Datasheet

Procalcitonin (PCT) Antibody HM365

Product Name	Procalcitonin (PCT) Antibody
Catalog Number	HM365
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	WB, LF, ELISA
Immunogen	Recombinant Procalcitonin
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgA
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.4 ± 0.2
Specificity	Procalcitonin
Matched Pairs	HM364 (detection)
Preservatives	0.05% Proclin 300