

Product Datasheet

Procalcitonin (PCT) Antibody HM365

| | |
|----------------------------|---|
| Product Name | Procalcitonin (PCT) Antibody |
| Catalog Number | HM365 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA, WB |
| Immunogen | Recombinant Procalcitonin |
| Brand | EastCoast Bio |
| Form/Appearance | Protein A Purified Monoclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype | IgA |
| Clonality | Monoclonal |
| Purity | >90% |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.4 |
| Specificity | Procalcitonin |
| Matched Pairs | HM364 (detection) |
| Preservatives | 0.05% Proclin 300 |