

Product Datasheet

Cystatin C Antibody HM367

Product Name	Cystatin C Antibody
Catalog Number	HM367
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	WB, LF, ELISA
Immunogen	Cystatin C native protein
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.4 ± 0.2
Specificity	Cystatin C
Matched Pairs	HM368 (detection)
Preservatives	0.05% Proclin 300