

Product Datasheet

Buprenorphine (BUP) Antibody HM376

Product Name	Buprenorphine (BUP) Antibody
Catalog Number	HM376
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	BUP-BSA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	> 90 %
Buffer	0.05M Phosphate, pH 7.15 with 0.85% Sodium Chloride
Specificity	BUP, norBUP and BUP glucuronide: 100%
Matched Pairs	LA283 or LA284 (capture)
Preservatives	0.05% Sodium Azide