

Product Datasheet

Buprenorphine (BUP) Antibody HM377

Product Name	Buprenorphine (BUP) Antibody
Catalog Number	HM377
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	BUP-BSA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	>90%
Buffer	0.05M Phosphate, pH 7.15 with 0.85% Sodium Chloride
Specificity	BUP, norBUP and BUP glucuronide: 100%
Matched Pairs	LA283 or LA284 (capture)
Preservatives	0.05% Sodium Azide