

Product Datasheet

Dehydroepiandrosterone sulfate (DHEA-S) Antibody HM380

Product Name	Dehydroepiandrosterone sulfate (DHEA-S) Antibody
Catalog Number	HM380
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	DHEA-BSA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2b, kappa
Clonality	Monoclonal
Purity	> 90 %
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl
Specificity	DHEA-3-sulfate 100%, DHEA 10%
Matched Pairs	LA285 (capture)
Preservatives	0.05% Sodium Azide