

Product Datasheet

EDDP Antibody HM381

Product Name	EDDP Antibody
Catalog Number	HM381
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	EDDP-KLH
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	> 90 %
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Specificity	EDDP 100%, Methadone 0.1%
Matched Pairs	LA282 (capture)
Preservatives	0.05% Sodium Azide