

Product Datasheet

Malaria (pLDH) Antibody HM382

Product Name	Malaria (pLDH) Antibody
Catalog Number	HM382
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Purified Malaria LDH
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	> 90 %
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Specificity	Malaria pLDH
Matched Pairs	HM383 (detection)
Preservatives	0.05% Sodium Azide