

Product Datasheet

Human TRAIL ELISA Kit 850.770.096

Product Name	Human TRAIL ELISA Kit
Catalog Number	850.770.096
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	93.75-3000 pg/ml
Sensitivity	64 pg/ml
Tested Applications	ELISA
Alternative Names	Apo2L, CD253
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human soluble molecules
Specificity	Recognizes both natural and recombinant soluble human Trail
Sample Types	Serum, Plasma, Cell culture supernatant