

Product Datasheet

Morphine Antibody HM385

Product Name	Morphine Antibody
Catalog Number	HM385
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Morphine-BSA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Meets 300 ng/ml cutoff
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	>90%
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Specificity	Morphine, Codeine, Hydrocodone 100%, Hydromorphone 90%
Matched Pairs	LA281 (capture)
Preservatives	0.05% Sodium Azide