

Product Datasheet

Pregnanediol Glucuronide (PDG) Antibody HM386

Product Name	Pregnanediol Glucuronide (PDG) Antibody
Catalog Number	HM386
Description	Purified Monoclonal Antibody
Tested Applications	ELISA, LF
Immunogen	PDG-BSA
Brand	EastCoast Bio
Form/Appearance	Purified monoclonal antibody
Storage	Short Term: 2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	≥ 90% by SDS-PAGE
Buffer	0.015 M potassium phosphate, pH 7.2 ± 0.2, with 0.85% sodium chloride, 0.05% sodium azide as a preservative.
Specificity	Pregnanediol glucuronide
Matched Pairs	LA286 (capture)