

Product Datasheet

Pregnanediol Glucuronide (PDG) Antibody HM387

Product Name	Pregnanediol Glucuronide (PDG) Antibody
Catalog Number	HM387
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	PDG-BSA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2b, kappa
Clonality	Monoclonal
Purity	> 90 %
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Specificity	PDG
Matched Pairs	LA286 (capture)
Preservatives	0.05% Sodium Azide