

Product Datasheet

Dengue Virus (DENV-NS1 Serotype 2,3,4) Antibody HM388

Product Name	Dengue Virus (DENV-NS1 Serotype 2,3,4) Antibody
Catalog Number	HM388
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Purified Dengue Virus
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	> 90 %
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Specificity	NS1 protein of Dengue Virus types 2,3,4
Matched Pairs	HM843 or HM349 (detection), HM389 (capture)
Preservatives	0.05% Sodium Azide