

Product Datasheet

Dengue Virus (DENV-NS1 Serotype 1) Antibody HM389

| | |
|----------------------------|--|
| Product Name | Dengue Virus (DENV-NS1 Serotype 1) Antibody |
| Catalog Number | HM389 |
| Description | Purified Monoclonal Antibody |
| Species/Host | Mouse |
| Tested Applications | ELISA, LF |
| Immunogen | Purified Dengue Virus |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Monoclonal Antibody |
| Storage | Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| Isotype | IgG2b, kappa |
| Clonality | Monoclonal |
| Purity | > 90 % |
| Buffer | 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride |
| Specificity | NS1 protein of Dengue Virus type 1 |
| Matched Pairs | HM349 (detection) and/or HM388 (capture) |
| Preservatives | 0.05% Sodium Azide |