

Product Datasheet

Ketamine Antibody HM390

Product Name Ketamine Antibody

Catalog Number HM390

Species/Host Mouse

Tested Applications LF, ELISA

Brand EastCoast Bio

Form/Appearance Purified Monoclonal Antibody

Storage Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

Isotype IgG2a, kappa

Clonality Monoclonal

Purity >90%

Buffer 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride

Cross Reactivity Ketamine 100%, Norketamine 60%

Specificity Ketamine

Matched Pairs LA280 (capture)

Preservatives 0.05% Sodium Azide