

Product Datasheet

Ketamine Antibody HM391

Product Name	Ketamine Antibody
Catalog Number	HM391
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2b, kappa
Clonality	Monoclonal
Purity	> 90 %
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	Ketamine 100%, Norketamine 60%
Specificity	Ketamine
Matched Pairs	LA280 (capture)
Preservatives	0.05% Sodium Azide