

## Product Datasheet

### Measles Virus Antibody HM403

|                            |  |
|----------------------------|--|
| <b>Product Name</b>        | Measles Virus Antibody   |
| <b>Catalog Number</b>      | HM403  |
| <b>Description</b>         | Purified Monoclonal Antibody   |
| <b>Species/Host</b>        | Mouse  |
| <b>Tested Applications</b> | ELISA, LF  |
| <b>Immunogen</b>           | Purified Measles Virus   |
| <b>Brand</b>               | EastCoast Bio  |
| <b>Form/Appearance</b>     | Purified Monoclonal Antibody   |
| <b>Storage</b>             | Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Monoclonal   |
| <b>Purity</b>              | > 95 % by SDS-PAGE   |
| <b>Buffer</b>              | 10 mM Phosphate Buffered Saline, pH 7.0                                      |
| <b>Specificity</b>         | Measles Virus nucleoprotein  |
| <b>Preservatives</b>       | 0.02% Sodium Azide   |