

Product Datasheet

Helicobacter pylori Antibody HM404

Product Name	Helicobacter pylori Antibody
Catalog Number	HM404
Description	Purified Monoclonal Antibody
Immunogen	Purified H. pylori
Brand	EastCoast Bio
Form/Appearance	Purified monoclonal antibody
Storage	Short Term: 2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	≥ 95 % by SDS-PAGE
Buffer	Phosphate-buffered saline, pH 7.0 ± 0.2, 0.09% sodium azide as a preservative.
Specificity	H. pylori
Matched Pairs	HM092 (detection) or HM405 (capture)