

# Product Datasheet

## Clostridium difficile toxin A Antibody HM409

|                     |                                                                           |
|---------------------|---------------------------------------------------------------------------|
| Product Name        | Clostridium difficile toxin A Antibody                                    |
| Catalog Number      | HM409                                                                     |
| Species/Host        | Mouse                                                                     |
| Tested Applications | LF, ELISA                                                                 |
| Immunogen           | Purified native C. difficile toxin A                                      |
| Brand               | EastCoast Bio                                                             |
| Form/Appearance     | Purified Monoclonal Antibody                                              |
| Storage             | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype             | IgG1                                                                      |
| Clonality           | Monoclonal                                                                |
| Purity              | >95% by SDS-PAGE                                                          |
| Buffer              | Phosphate Buffered Saline, pH 7.0                                         |
| Specificity         | C. difficile toxin A                                                      |
| Matched Pairs       | HM408 (capture)                                                           |
| Preservatives       | 0.09% Sodium Azide                                                        |